

Best in the Microworld, Better beyond the Globe



ValvuloPro™

Valvuloplasty Balloon Catheter

**Valvulo Pre-dilation & Delivery
Strength & Precision Expansion
Anti-slip Structure**



TT Medical, Inc.
www.ttmedicalinc.com



TT Medical, Inc.

ValvuloPro™ Valvuloplasty Balloon Catheter Features

ValvuloPro™ Balloon Dilator is intended to be used to pre-dilate heart valves including aorta, mitral and tricuspid valves and delivery the transcatheter heart valves.

★ **Precision Dilatation:** Precision controlled balloon dimension during expansion.

★ **Strength Dilatation:** The balloon was made with optimized high pressure material for strength and rupture resistant.

★ **Anti-slip** Design enhanced balloon dilation performance and valve deliverability.

★ **Quick Inflation and Deflation:** the optimized dual lumen catheter design make the balloon inflation and deflation quick.



www.ttmedicalinc.com



TT Medical, Inc.

ValvuloPro™

Valvuloplasty Balloon Catheter

Anti-slip Design

Uniform Expansion

Strength & Precision
Dilation

Quick Inflation & Deflation



www.ttmedicalinc.com



Order Information (For OEM only) :

Number Reference	Inflated Balloon Diameter (mm)	Balloon Length (mm)	Dilatation Pressure (ATM)	RBP (ATM)
VBC1540	15	40	3	6
VBC1840	18	40	3	6
VBC2140	21	40	3	6
VBC2340	23	40	3	6
VBC2640	26	40	3	6
VBC2940	29	40	3	6

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

ValvuloPro™ is a trademark or registered trademark of TT Medical. All cited trademarks are the property of their respective owners.

